#### Conceptual Chemistry Introduction to Molecular Spectroscopy: Glow In The Dark Explorations with Tonic Water and Colorimetric Analysis of Tabletop Sweeteners

Christopher J. FenkHeathClaudia Khourey-BowersDonalDerek FaielloDan BCarol ThombsGretchGene EasterGretch

Heather Aubihl Donald G. Gerbig, Jr. Dan Brubaker Gretchen Kovack-Pappas

Kent State University at Tuscarawas/Stark

### **Conceptual Chemistry**

**Conceptual Chemistry** is a graduate course designed for grade school and middle school teachers to assist in their understanding of chemistry and to provide concrete ideas that they can take back to their classrooms to teach their students.

# **Conceptual Chemistry**

Participants in this course receive:

- Free tuition and five graduate credit hours from the College of Education of Kent State University. (\$2,340 value)
- Over \$850 worth of materials and supplies to take back to the classroom.

# **Conceptual Chemistry**

Support for *Conceptual Chemistry* and the development/production of this material was provided by a grant under the federally funded *Improving Teacher Quality State Grants Program*, administered by the Ohio Board of Regents.

# Spectroscopy

#### Spectroscopy

The investigation into the nature of matter using electromagnetic radiation.

### Most Common Types

Infrared, ultraviolet, visible, NMR, X-ray and microwave.

# Activity Objectives

### Key concepts:

- Molecular Spectroscopy
- Chemical Changes
- Physical Changes
- Quantitative and Qualitative Analyses
- Having fun with science!

#### Purpose

To determine the concentration of drink mix in an unknown sample using semi-quantitative visible light spectroscopy.





### Spectroscopic Analysis

#### Procedure

- 1. Obtain a 15 ml plastic conical tube containing Crystal Light drink mix (*12 ml "stock" solution*).
- 2. Use a plastic pipette to withdraw 8.0 ml of Crystal Light drink mix from the conical tube containing the "stock" solution and place it into the <u>first</u> test tube. (*This leaves 4.0 ml drink mix in the plastic tube.*)

### Spectroscopic Analysis

#### Procedure (cont.)

- 3. Add 8.0 ml water to the plastic conical tube, replace the cap and mix thoroughly. (*This gives 12 ml of solution again!*)
- 4. Withdraw 8.0 ml of the diluted Crystal Light drink mix from the conical tube and place it into the <u>second</u> test tube. (*This leaves 4.0 ml drink mix in the plastic tube again.*)

### Spectroscopic Analysis

#### Procedure (cont.)

- 5. Add 8.0 ml water to the plastic conical tube, replace the cap and mix thoroughly.
- Continue this dilution process three more times until <u>five</u> glass test tubes are filled.



### Spectroscopic Analysis

### Procedure (cont.)

- Obtain a conical tube containing the *"unknown"* sample. Place 8.0 ml of the *"unknown"* into the <u>sixth</u> test tube.
- 8. Compare the *"unknown"* sample to the five standards to determine the dilution ratio of the unknown sample. (1.9 g drink mix dissolved into 500 ml water).



|   | Spectroscopic Analysis |         |          |          |           |  |  |  |  |
|---|------------------------|---------|----------|----------|-----------|--|--|--|--|
| Dilution Math   |                        |         |          |          |           |  |  |  |  |
| $\frac{3.8 \text{ g/}}{\text{L}} \times \frac{4 \text{ ml}}{12 \text{ ml}} = \frac{1.3 \text{ g/}}{\text{L}}$ |                        |         |          |          |           |  |  |  |  |
|   | Tube #1                | Tube #2 | Tube #3  | Tube #4  | Tube #5   |  |  |  |  |
|   | 3.8 g/L                | 1.3 g/L | 0.43 g/L | 0.14 g/L | 0.058 g/L |  |  |  |  |
| "Unknown" was prepared to match "Tube #2"   |                        |         |          |          |           |  |  |  |  |



#### Materials

12 ml tonic waterDistilled water6 glass test tubes1.0 ml graduate plastic pipetUV light source

Tonic water is used as purchased and assumed to contain 60 mg/L quinine

### Spectroscopic Analysis

### Procedure

- 1. Obtain a 15 ml plastic conical tube containing tonic water (*12 ml "stock" solution*).
- 2. Use a plastic pipette to withdraw 8.0 ml of tonic water from the conical tube containing the "stock" solution and place it into the <u>first</u> test tube. (*This leaves 4.0 ml tonic water solution in the plastic tube.*)

### Spectroscopic Analysis

### Procedure (cont.)

- 3. Add 8.0 ml water to the plastic conical tube, replace the cap and mix thoroughly. (*This gives 12 ml of solution again!*)
- Withdraw 8.0 ml of the diluted tonic water from the conical tube and place it into the <u>second</u> test tube. (*This leaves 4.0 ml drink mix in the plastic tube again.*)

#### Procedure (cont.)

- 5. Add 8.0 ml water to the plastic conical tube, replace the cap and mix thoroughly.
- Continue this dilution process three more times until <u>five</u> glass test tubes are filled.
- Shine a UV light onto the test tubes being careful not to expose yourself or others to the light. Note the variable brightness of the tubes with respect to concentration.

### Spectroscopic Analysis

#### Procedure (cont.)

- Obtain a conical tube containing the *"unknown"* sample. Place 8.0 ml of the *"unknown"* into the <u>sixth</u> test tube.
- Compare the "unknown" sample to the five standards to determine the dilution ratio of the unknown sample. (60 ppm quinine is dissolved into 1,000 ml solution).

### Spectroscopic Analysis

#### **Dilution Math**

Tonic water is regulated and may only contain up to 83 ppm of quinine. Most commercial tonic water ranges from 25 -60 ppm.

$$60 \text{ ppm x } \frac{4 \text{ ml}}{12 \text{ ml}} = 20 \text{ ppm}$$

| Tube #1 | Tube #2 | Tube #3 | Tube #4 | Tube #5  |
|---------|---------|---------|---------|----------|
| 60 ppm  | 20 ppm  | 6.7 ppm | 2.2 ppm | 0.70 ppm |

"Unknown" was prepared to match "Tube #3"

## History of Tonic Water

Quinine is an anti-pyretic, anti-inflammatory, analgesic, anti-malarial bitter tasting drug and has been used for hundreds of years to treat malaria symptoms.





## Ultraviolet Light

### UV-A

320-400 nm, lowest energy UV light, suspected of damaging skin (premature aging and wrinkling), and skin cancer.



### **Modified Biuret Reagent**

Biuret Reagent (250 mL): Copper (II) sulfate pentahydrate (3.75 g, 15.0 mmol) was added to a 250 mL volumetric flask. Water was added (40 mL) and the copper (II) sulfate pentahydrate was completely dissolved. A solution of potassium sodium tartrate tetrahydrate (18.82 g, 52.5 mmol) in water (60 mL) was prepared in a separate container and combined with the copper (II) sulfate solution. Finally, dilute the mixture to volume using 1  $\underline{M}$  sodium carbonate solution.

## Modified Biuret Reagent

Sodium Carbonate (1 M): Sodium carbonate (105.99 g, 1 mol) was added to a 1 L volumetric flask and dissolved into 750 mL water. The solution was then diluted to volume with water.



# **Qualitative Analysis**

1. Use dropper bottle to add 20 drops of reference materials to 24 well plate.

Addition Order

- B1: water, B2: aspartame solution
- 2. Add 5 drops biuret reagent.
- 3. Note results.









Splenda® – sucralose, dextrose and maltodextrin













For 14 trials:

For 44 Equal® Samples:

aspartame per sample = 37.6 ± 2 mg/g

manufacturer' s claim = 38 mg/g

1.05% difference!

# Conclusions

- New method for identifying and determining the amount of aspartame in tabletop sweeteners.
- Analysis of 44 samples gave 37.6 mg ± 2 mg aspartame/g Equal<sup>®</sup> observed compared to 38 mg/g as stated by the manufacturer. (1.05% difference)
- Christopher J. Fenk, Donald G. Gerbig, Jr. and Nathan Kaufman "A New Colorimetric Assay of Tabletop Sweeteners Using a Modified Biuret Reagent" J. Chem. Educ. 2007, 84 (10), 1676-1678.

## Conclusions

- Ø Methods for introducing the concepts of spectroscopic analysis of samples were presented.
- Ø Both visible and UV light were used to semiquantitatively determine the concentration of unknown samples.
- Ø Unknowns are determined using the oldest/ most sophisticated detection device known. Your eyes.